

An English guide to experiments



Heavy Ion Nuclear Physics Laboratory

GIBELIN Julien

March 29, 2010

Forewords

This simple help has been created as an introduction for analysis at RIKEN CYCLOTRON and especially RIPS (E6 room), the RIKEN Projectile fragment Separator. The author hopes than incomers in RIKEN CYCLOTRON will find help basis during the course of experiments.

Writing conventions are :

< > compulsory parameter

{[]} optional parameter

/ choice

-> ‘ ‘ ’ ’ keyboard input

Examples are signaled by $\diamond\diamond$.

Please send all comments, bugs... to the author :

gibelin@rarfaxp.riken.go.jp

A *pdf* version of this booklet is available for downloading at:

http://rarfaxp.riken.go.jp/~gibelin/riken_guide.pdf

A *html* version of this booklet is also available at:

http://rarfaxp.riken.go.jp/~gibelin/riken_guide/

Last version: March 29, 2010

Acknowledgments

Especially to the author of ANAPAW & the data acquisition system used by the Heavy Ion Laboratory (Heavion): S. TAKEUCHI & H. BABA.

My others colleagues played a great role in this guide by correcting my numerous English mistakes !

Merci !

J. GIBELIN

PS: Almost 5 years later I am back to RIKEN for an experiment, so I took this opportunity to update this manual, with the help of my colleagues from LPC Caen.

Contents

1	Data Acquisition System	1
1.1	Basis of DAQ system in RIKEN	1
1.2	Bbcaenhv	2
1.2.1	Header overview	2
1.2.2	Command description	2
	Crate	2
	Edit	2
	Load	2
	Save	2
	On/Off	2
	Update	2
	Monitor	3
	Quit	3
1.3	BabarDAQ program's commands	3
1.3.1	MOUNT	3
1.3.2	MTON	3
1.3.3	HDON	3
1.3.4	WTH	3
1.3.5	START	3
1.3.6	NSSTA	3
1.3.7	STOP	3
1.3.8	PRISCA	4
1.3.9	MTOFF	4
1.3.10	DISMOU	4
1.3.11	PRIHEADER	4
1.3.12	SETNUMBER	4
1.3.13	Main procedure	4
1.4	BabarDAQ Emergency Stop	5
	To restart	5
1.5	SSDHV	5
1.5.1	GETFILE	6
1.5.2	PUTFILE	6
1.5.3	SET	6
1.5.4	GET	6
1.5.5	CTRL	6
1.5.6	ZERO	6
1.5.7	NAME	6
1.5.8	MONITOR	6

CONTENTS

1.5.9	PRINT	6
1.5.10	EXIT	7
1.5.11	STAT	7
1.5.12	VIEW	7
1.5.13	HELP	7
1.5.14	Unix valuable commands	7
2	Anapaw	9
2.1	Install ANAPAW	9
2.1.1	ANAPAW versions	9
2.1.2	Before installing	9
	Off line analysis	10
	On line analysis	10
2.1.3	How to log in user directory	10
	Off line analysis	10
	On line analysis	11
2.1.4	Compiling ANAPAW	11
2.2	Software	11
2.2.1	Login	11
2.2.2	ANAPAW structure	12
2.3	Commands	12
2.3.1	/ANALYS/	12
	/ANALYS/LOOP	12
	/ANALYS/STATUS	13
	/ANALYS/RDMP	13
	/ANALYS/ANTUPLE	14
2.3.2	/ANALYS/DEFINITION	14
	/ANALYS/DEFINITION/BOOK	14
	/ANALYS/DEFINITION/CLEAR	14
	/ANALYS/DEFINITION/ANAADD	14
	/ANALYS/DEFINITION/AWRITE	15
	/ANALYS/DEFINITION/CHGATE	15
	/ANALYS/DEFINITION/CHHST	16
	/ANALYS/DEFINITION/HST1	16
	/ANALYS/DEFINITION/HST2	17
	/ANALYS/DEFINITION/LIST1D	17
	/ANALYS/DEFINITION/LIST2D	18
	/ANALYS/DEFINITION/LISTPROF	18
	/ANALYS/DEFINITION/LISTGATE	18
	/ANALYS/DEFINITION/AVIEW	18
	/ANALYS/DEFINITION/HCUT	19
2.3.3	/ANALYS/HISTOGRAM/	19
	/ANALYS/HISTO/HSTORE	19
	/ANALYS/HISTO/FETCH	19
	ANALYS/HISTOGRAM/HLIST or I	20
	/ANALYS/HISTO/HT	20
	/ANALYS/HISTO/HTP	20
	/ANALYS/HISTO/HNT or /ANALYS/HISTO/HN	20
	/ANALYS/HISTO/HBT or /ANALYS/HISTO/HB	20
	/ANALYS/HISTO/BLOW	21

CONTENTS

	/ANALYS/HISTO/XBLOW	21
	/ANALYS/HISTO/XYBLOW	21
	/ANALYS/HISTO/PRX	21
	/ANALYS/HISTO/PRY	21
	/ANALYS/HISTO/BNX	21
	/ANALYS/HISTO/BNY	21
	/ANALYS/HISTO/SLX	22
	/ANALYS/HISTO/SLY	22
	/ANALYS/HISTO/MAMI	22
	/ANALYS/HISTO/XSTATUS	22
	/ANALYS/HISTO/HSTATUS	22
	/ANALYS/HISTO/ERASE	22
	/ANALYS/HISTO/HDELETE	22
	/ANALYS/HISTO/XVAL; /ANALYS/HISTO/XYVAL	23
	/ANALYS/HISTO/HSTAT	23
	/ANALYS/HISTOGRAM/XSTAT	23
	/ANALYS/HISTOGRAM/FIGA	23
	/ANALYS/HISTOGRAM/CPRO	23
	/ANALYS/HISTOGRAM/AVY	23
	/ANALYS/HISTOGRAM/XFITG	24
2.3.4	/GRAPHICS/	24
	/GRAPHICS/PSON	24
	/GRAPHICS/PSOFF	24
3	Analys Code	25
3.1	Hierarchy	25
3.2	Comments	25
3.3	Analys	26
3.4	Gates	26
3.4.1	GATE	26
3.4.2	Logical OR, AND	26
3.4.3	XYGATE	26
3.4.4	STOP	27
3.5	Histograms	27
3.5.1	HST1: 1D Histograms	27
3.5.2	PROFILE	27
3.5.3	HST2: 2D Histograms	28
3.6	ANA-CODE file example	28
4	Anapaw Source Code	31
4.1	Analyzer Definition	31
4.2	Encoding File	32
4.2.1	From raw data to observable	32
4.2.2	A first example	32
4.2.3	Calibration	33
4.2.4	Multi array detector	35
4.3	Ntuples	36

CONTENTS

5	Miscellaneous	39
5.1	Phone numbers	39
5.2	Data copy	39
5.2.1	Tape duplication	39
5.2.2	From tape to hard drive	39
5.2.3	From hard drive to hard drive	39
5.3	Data structure	40
5.4	Scaler data	41

Chapter 1

Data Acquisition System

Data Acquisition System system for Heavy Ion Nuclear Physics Laboratory has been mainly supervised by H. BABA. For more details refer to his home-page (in Japanese): <http://daq.rikkyo.ac.jp/babar/>

1.1 Basis of DAQ system in Riken

For DAQ, both WME and CAMAC are used. Below is a simplified scheme¹.




Figure 1.1: Typical configuration of BabarDAQ. Some client PCs can lie on the Ethernet

Here, our concern is not on how data are read from electronic module but mainly on how we can *monitor* them. So as let's describe three main programs : BBCEAHV (see Section 1.2) which allows you to read and set high voltage

¹Extract from *Development of New Data Acquisition System BabarDAQ*, RIKEN Accel. Prog. Rep. 34 (2001)

1.2. Bbcaenhv

of some detectors, BABARL which allows you to launch DAQ and SSDHV for monitoring silicons H.V.

First, let's give name to current main computers² :

rips1 controller

ripscc CAMAC control

rips00 analysis

bbhv H.V. control

Finally, you can find in the miscellaneous part (see Section 5.3), the complete description of data structure.

1.2 Bbcaenhv

Allow you to read and edit high voltage power supply values. You must first go to the H.V. computer control (mainly bbhv, ask for password).

1.2.1 Header overview

CRATE ## PAGE #/# Crate Page Edit Load Save On/Off Update Monitor Quit

1.2.2 Command description

Crate

Change/read crate number. Keyboard command : **C**

Edit

Edit selected field. Keyboard command : **E**

Load

Load H.V. configuration's file. Keyboard command : **L**

Save

Save H.V. configuration's file. Keyboard command : **S**

On/Off

Turn **On/Off** the HV for selected device. Keyboard command : **O**

Update

Update values. Keyboard command : **U**

²Note : these appellations might alter from time to time...

1.3. BABARLDAQ PROGRAM'S COMMANDS

Monitor

Update values every 10 seconds. Keyboard command : **M**

Quit

Exit BBCEAHV. Keyboard command : **Q**

1.3 BabarLDAQ program's commands

1.3.1 MOUNT

Mount physically the tape on hardware.

```
Command> MOUNT
```

1.3.2 MTON

Mount the **t**ape **o**n software.

```
Command> MTON
```

1.3.3 HDON

Use **h**ard **d**rive to write **o**n.

```
Command> HDON
```

1.3.4 WTH

Write the **h**eaders of the next run.

```
Command> WTH
```

1.3.5 START

Start a saved run.

```
Command> START
```

1.3.6 NSSTA

Non saved **s**tart.

```
Command> NSSTA
```

1.3.7 STOP

Stop the current run and ask for lender.

```
Command> STOP
```

1.3. BABARLDAQ PROGRAM'S COMMANDS

1.3.8 PRISCA

Print the CAMAC **scaler** values.

```
Command> PRISCA
```

1.3.9 MTOFF

Unmount the tape **off** from software, swapping automatically to hard drive (see Section 1.3.3).

```
Command> MTOFF
```

1.3.10 DISMOU

Unmount/eject the tape from hardware.

```
Command> DISMOU
```

1.3.11 PRIHEADER

Print the **header** and the **ender** of a raw data file.

```
Command> PRIHEADER <FNAME>
```

1.3.12 SETNUMBER

Change run number (next RUN is *set number + 1*).

```
Command> SETNUMBER <NUM>
```

1.3.13 Main procedure

◆◆ If you want to make a non recorded acquisition :

```
Command> nssta
```

◆◆ If you want to make a recorded acquisition :

```
Command> mount
```

Mount physically the tape on hardware.

```
Command> mton
```

Mount the tape, case of software

```
Command> wth
```

```
Header :
```

```
-> ‘ ‘ Bonjour ceci est une acquisition ’ ’
```

```
Command> start
```

... And wait for you data ...

```
Command> stop
```

```
Ender :
```

```
-> ‘ ‘ Ceci etait une acquisition ’ ’
```

1.4. BABARLDAQ EMERGENCY STOP

Stop the acquisition and write the ender

```
Command> prisca
```

...and don't forget to print the scaler values ! After that you can restart another acquisition. If you want to change the tape.

```
Command> mtoff
```

```
Command> dismou
```

1.4 BabarLDAQ Emergency Stop

- End the program using **Ctrl+C**
- In *rips1* shell, list process using **ps x**.
- ...
- Change the console by using **Alt+F2**
- In
- *ripscc* shell, run **ps ax**. Beware that this command will list all processes in the terminal and is thus non appropriated on *rips1* You should then see **bbcamac**, **collector**, and **transfer**. Kill only **transfer** while logged as root.
- Run **/sbin/lsmmod**. You should see a module named **babarldaq**. Kill it by running **/sbin/rmmod babarldaq**, as well as all other **nddq*** or **bb-*** modules.
- ...
- Go back to the first console: **Alt+F1**
- Run **netstat -a | less**. Look for any '167.***' tcp number. ?? If processes remain **TIME_WAIT** for more than 1 minute, this means that they (and especially **transfert**) remain in *ripscc* or *rips1*, so please go back and kill them.
- ??

To restart

Log to

SSDHV allows you to monitor and control HV crate for Silicon detectors. Before lunching SSDHV program you should ssh to **gomigo@gomigo** (ask for the password) and then lunch SSDHV.

1.5 SSDHV

Monitor for silicon H.V. Lunch with **ssdhv**, login **ssdhv**.

1.5. SSDHV

1.5.1 GETFILE

Get a status file

```
SSDHV : GETFILE <FNAME>
```

1.5.2 PUTFILE

Put status into file

```
SSDHV : PUTFILE <FNAME>
```

1.5.3 SET

Set voltage (in V), current (in nA) for a specified channel, increasing the voltage with a specified speed (V/s).

```
SSDHV : SET <CHANNEL> <VOLTAGE> <CURRENT> <SPEED>
```

1.5.4 GET

Get/update information on the SSD status.

```
SSDHV : GET [CHANNEL]
```

1.5.5 CTRL

```
SSDHV : CTRL <ON/OFF>
```

1.5.6 ZERO

Set the high voltage to zero for *all* channel. Slowly of course !.

```
SSDHV : ZERO
```

1.5.7 NAME

Give a channel a name. Only visible after updating (see Section [1.5.4](#))

```
SSDHV : NAME <CHANNEL> [CHNAME]
```

Warning: no possible space, example :

```
SSDHV : name 1 ssd_1-1
```

1.5.8 MONITOR

Update values every 10 seconds.

```
SSDHV : MONITOR <ON/OFF>
```

1.5.9 PRINT

Print SSD status.

```
SSDHV : PRINT
```

1.5.10 EXIT

Exit.

```
SSDHV : EXIT
```

1.5.11 STAT

Print out statistics on a file.

```
SSDHV : STAT <FNAME>
```

1.5.12 VIEW

View the contents of a file.

```
SSDHV : VIEW <FNAME>
```

1.5.13 HELP

View help.

```
SSDHV : HELP
```

1.5.14 Unix valuable commands

You can use LS, CD & PWD.

Chapter 2

Anapaw

ANAPAW , made by S. TAKEUCHI, combines analysis (in real time) with PAW. For further informations, please consult his site (in japanese) :

<http://www.ne.rikkyo.ac.jp/~takeuchi/>

2.1 Install Anapaw

2.1.1 Anapaw versions

Because ANAPAW uses CERN Paw and CERNLIB, you should first correctly install this last softwares, see

<http://it-div.web.cern.ch/it-div/physicscomputing/>

Two versions exist : the latest and fastest is

<http://rarfaxp.riken.go.jp/~takesato/anapaw/Download/download.html> whereas the oldest one is

<http://rarfaxp.riken.go.jp/~takesato/anapaw/Download/download.html>, the main difference lies in the way the software deals with memory.

The analysis part is possible thanks to software named ANALYS. Basic source codes should be download for it. The latest version uses :

<http://rarfaxp.riken.go.jp/~takesato/anapaw/Download/download.html>, the other :

<http://rarfaxp.riken.go.jp/~takesato/anapaw/Download/download.html>

2.1.2 Before installing

Warning: ANAPAW commands have been written using C SHELL, thus you should swap to this shell.

Edit or create the `.cshrc` file at the root of your home directory.

If we suppose that the ANAPAW source files are in `~/usr/anapaw` and you want your work directory to be `~/exp/r300n`, add the following:

```
#ANAPAW
alias analogin 'source $HOME/usr/anapaw/Setup/setupanapaw/'
setenv ANAPAW_WORK $HOME/exp/r300n
```

Then for this change to be effective :

2.1. INSTALL ANAPAW

```
source .cshrc
```

ANAPAW uses a set of environment variables, contained in the `setupanapaw` file. At least two of them must be changed to fit to your computer's configuration. First, as far as ANAPAW deals CERN's PAW - and other libraries - specify the location of CERNLIB. In second the ANAPAW source files location itself should be specified.

As an example, the first lines of `setupanapaw` can look like (defaults):

```
setenv CERNLIB      /cern/pro/lib
setenv ANAPAW_HOME $HOME/usr/local/anapaw
```

Off line analysis

ANAPAW can be use for on-line analysis as well as off line analysis.

The suitable configuration for local and off-line is :

```
# For using on local.
source $ANAPAW_SETUP/setup.local
```

You must create work directory, for examples, if you want to to work in `~/exp/r300n`, you must create : `~/exp`. Or if you want `~/exp/r300n/offline` then create `~/exp/r300n`

Finally copy source files (basic ones are in `usersrc-*.tar.gz`) to e.g `~/exp/r300n/src` (or `~/exp/r300n/offline/src`). This can be off course change in the `setupanapaw` file (see Section 2.1.2).

On line analysis

For the on-line case, more than one user can make analysis, but everyone must use, at the beginning, the same source code. Then if the main directory will be `~/exp/r300n` :

```
> mkdir ~/exp/
> mkdir ~/exp/r300n
> mkdir ~/exp/r300n/src
> mkdir ~/exp/r300n/users
```

The `analogin` (see Section 2.1.3) commands will automatically create users files.

2.1.3 How to log in user directory

Off line analysis

After definition of main parameters in `setupanapaw` (see Section 2.1.2) you simply have to typewrite :

```
> analogin
```


On line analysis

Due to several user, you will have to specified your login. $\diamond\diamond$ For example for the user named *hogehoge*.

```
> analogin

== ANAPAW_SETUP ==

Please Input Your Name.

-> 'hogehoge'
```

Then if non directory exists for you, `analogin` will create it.

2.1.4 Compiling Anapaw

You first have to build the library, independently of the characteristics of the experiment.

Supposing that the `ANAPAW_SOURCE` variable has been assigned by the 'analogin' program (see Section 2.1.2). Just type :

```
> make -C $ANAPAW_SOURCE clean
> make -C $ANAPAW_SOURCE
> make -C $ANAPAW_SOURCE install
```

This series of commands can be simplified by using

```
> makelib
```

Then you must build the specific encoding files for the analysis. After changing to the directory where the user source files are (example : `~/exp/r300n/src`), type :

```
> make clean
> make
```

Everything can be simplified by using (from anywhere)

```
> makeana
```

2.2 Software

2.2.1 Login

Simply typewrite `analogin` and follow the instruction. Most of time, except when you want to create a personal ANAPAW setup, your login is `rips`.

Note that you can run a kumac in batch using `-b` argument, followed by the macro's name

```
rips@host > analogin
```

2.3. COMMANDS

2.2.2 Anapaw structure

The purpose here is not to explain in detail the ANAPAW structure, but to give a rough idea.

2.3 Commands

2.3.1 /ANALYS/

/ANALYS/LOOP

Event loop shell : this command allows you to read raw data, so as to analyze them. You can do it OFF-LINE or ON-LINE.

```
ANAPAW> ANA/LOOP [ FNAME CHOPT ]

      FNAME : ' ' On line reading
              Analys filename
      CHOPT  : ' ' Command mode
              B  Batch mode
```

◇◇Let's start with an example of simple OFF-LINE analysis, using the raw data contained in **run1018.rdf** and the analysis command of **sample.ana**

```
ANAPAW> book sample.ana
```

Read the analysis command from **sample.ana**, see Section [2.3.2](#)

```
ANAPAW> hst1
HST1> (=0)
```

1D-Histograms:

```
 1 ( HID : 101 ) ( TITLE : ch1                ) ( GATE : 0 )
   ( bin : 200 ) ( xmin : 0.00 ) ( xmax : 2000.00 )

 2 ( HID : 102 ) ( TITLE : ch2                ) ( GATE : 0 )
   ( bin : 200 ) ( xmin : 0.00 ) ( xmax : 2000.00 )

 3 ( HID : 103 ) ( TITLE : ch3                ) ( GATE : 0 )
   ( bin : 200 ) ( xmin : 0.00 ) ( xmax : 2000.00 )

 4 ( HID : 104 ) ( TITLE : ch4                ) ( GATE : 0 )
   ( bin : 200 ) ( xmin : 0.00 ) ( xmax : 2000.00 )
```

Default 1D-histogram proposed, see Section [2.3.2](#)

```
ANAPAW> ana/loop run1018.rdf
```

Open the data file **run1018.rdf**

```
ANAPAW/EVTLOOP>
ANAPAW/EVTLOOP> start
```

2.3. COMMANDS

Start reading the values.

```
ANAPAW-M : Interrupt Event Loop!  
ANAPAW/EVTLOOP>
```

This message appears when hitting , so as to stop the reading.

```
ANAPAW/EVTLOOP> scat  
SCAT> on  
ANAPAW/EVTLOOP>
```

This allows you to open a scatter histogram window. Especially useful to view on time 2D histograms.

```
ANAPAW/EVTLOOP> quit
```

Quit !

◇◇How to make ONLINE analysis ? Simply write:

```
ANAPAW> ana/loop
```

/ANALYS/STATUS

```
ANAPAW> ana/status  
  
ANAPAW-M : Current Status.  
  
Online Mode  
  
Blocks      : 0  
Valid Events : 0  
  
RDF :  
ANA : ana/gcheck.ana  
  
Analyzer : 1  
Analyzer : 3  
  
Now on Main Prompt.
```

/ANALYS/RDMP

Store the raw data to a file with a given gate. For gate help see Section 3.4. Note that the file has to be opened before any kind of storage !

```
ANAPAW> rdmp [ OC filename GID ]  
  
OC      : 0 Open file  
        C Close file  
filename : Raw data file name  
GID      : Gate IDentifier
```

2.3. COMMANDS

◇◇Example :

```
ANAPAW> rdmp 0 test.rdf
ANAPAW> ana/loop
```

And after taking data ...

```
ANAPAW> rdmp C
```

/ANALYS/ANTUPLE

Store the raw data to N-tuple files with a given gate.

```
ANAPAW> antuple [ OC GID RunNum ]

OC          :  O Open file
             :  C Close file
GID         :  Gate IDentifier
RunNum      :  Run number
```

Warning: Please note that this function needs special haking in `add_ntuple.f` file, see Section 4.3 for more details.

2.3.2 /ANALYS/DEFINITION

/ANALYS/DEFINITION/BOOK

Open analysis input file and book histograms. This procedure must be call before taking or reading data.

```
ANAPAW> book [ FNAME ]
```

FNAME : File name

/ANALYS/DEFINITION/CLEAR

Clear all histogram definitions, histograms and event summary.

```
ANAPAW> clear
```

/ANALYS/DEFINITION/ANAADD

Add analysis parameter.

```
ANAPAW> anaadd [ filename ]
```

◇◇Here under are simple commands for anaadd

```
ANAPAW> anaadd test.ana
FileName :
test.ana
```

Type 'AQ' to quit ANAADD mode.

```
Analys Com >
Analys Com >del
```

Delete last line.

```
Analys Com >aq
```

Exit from editor, save file and update.

/ANALYS/DEFINITION/AWRITE

Write analysis parameters to a given filename.

```
ANAPAW> awrite [ filename ]
```

◇◇Example:

```
ANAPAW> awrite anacode.ana
```

/ANALYS/DEFINITION/CHGATE

Change definition of given gate, as describe in ANA-CODE. You can change definition of [GATE], [AND] and [OR].

```
ANAPAW> chgate [ GID ]
```

```
      GID      : Gate IDentifier in ANA-CODE
```

◇◇◇Example of a simple GATE and a gate using 2 gates logically linked by OR.

```
ANAPAW> chgate 1
```

```
Kind : GATE
```

```
      ( GID :    1 ) ( KIND : Gate )
      ( VALUE :    1    1    1    4 ) ( Limit :    0.50    1.50 )
      ( TOTAL ENTRIES :                0 ) ( ACCEPTED ENTRIES :                0 )
```

```
Analyzer > (=1)
```

```
Word1    > (=1)
```

```
ANAPAW> chgate 11
```

```
Kind : AND
```

```
      ( GID :   11 ) ( KIND : And )
      ( ACCEPTED EVENTS :                0 )
```

```
ELEMENTS :
```

```
9 10
```

```
Elements > (<CR>=No Change)
```

2.3. COMMANDS

/ANALYS/DEFINITION/CHHST

Change definition of a given histogram, as described in ANA-CODE, see Section 3. For 1D histograms see Section 3.5.1, and for 2D histograms see Section 3.5.3.

```
ANAPAW> chhst [ HID CHOPT ]
```

```
HID   : Histogram ID to change a definition.
CHOPT : Options
      ' ' = Change a definition by Command Line.
      X   = Change a X-Range on X-window by mouse.
      Y   = Change a X-Range on Y-window by mouse.
      T   = Change a XY-Range on X-window by mouse.
```

◇◇Example:

```
ANAPAW> chhst 287
```

```
1D-Histograms:
```

```
( HID : 287 ) ( TITLE : Multiplicity           ) ( GATE :      0 )
( bin :  11 ) ( xmin :   -0.50 ) ( xmax :   10.50 )
( VALUE :   4 200 200  2 )
```

```
GateID  > (<CR>=0)
Analyzer > (<CR>=4)
ID 1    > (<CR>=200)
ID 2    > (<CR>=200)
Word    > (<CR>=2)
Bin     > (<CR>=11)
Min     > (<CR>=-0.5)
Max     > (<CR>=10.5)
Title   > (<CR>=Multiplicity)
```

```
ANAPAW> chhst 1 x
```

/ANALYS/DEFINITION/HST1

Define ZONE ,IDs of 1D-histograms to see. By defaults (press or value given) the output is created with the first 4 histograms disposed 2 × 2

```
ANAPAW> hst1 [ NX NY ]
```

```
NX : Number of divisions along X
NY : Number of divisions along Y
```

◇◇Example of default result :

```
ANAPAW> hst1
HST1> (=0)
```

```
1D-Histograms:
```

2.3. COMMANDS

```
( HID : 101 ) ( TITLE : ch1 ) ( GATE : 2 )
( bin : 2 ) ( xmin : -0.50 ) ( xmax : 1.50 )
( VALUE : 1 1 1 5 )

( HID : 102 ) ( TITLE : ch2 ) ( GATE : -2 )
( bin : 2 ) ( xmin : -0.50 ) ( xmax : 1.50 )
( VALUE : 1 1 1 5 )

( HID : 103 ) ( TITLE : ch3 ) ( GATE : 0 )
( bin : 11 ) ( xmin : -0.50 ) ( xmax : 10.50 )
( VALUE : 3 100 100 2 )

( HID : 104 ) ( TITLE : ch4 ) ( GATE : 3 )
( bin : 11 ) ( xmin : -0.50 ) ( xmax : 10.50 )
( VALUE : 3 100 100 2 )
```

◇◇Example 2, where the defaults system allows you to simply replace the first histogram.

```
ANAPAW> hst1
HST1> (=0) 105
HST1> (=105) 0
```

1D-Histograms:

```
( HID : 105 ) ( TITLE : Coin 2 Gated ) ( GATE : 2 )
( bin : 2 ) ( xmin : -0.50 ) ( xmax : 1.50 )
( VALUE : 1 1 1 5 )

( HID : 106 ) ( TITLE : Coin 2 Gated inv ) ( GATE : -2 )
( bin : 2 ) ( xmin : -0.50 ) ( xmax : 1.50 )
( VALUE : 1 1 1 5 )

( HID : 107 ) ( TITLE : M_dE ) ( GATE : 0 )
( bin : 11 ) ( xmin : -0.50 ) ( xmax : 10.50 )
( VALUE : 3 100 100 2 )

( HID : 108 ) ( TITLE : M_dE Gated ) ( GATE : 3 )
( bin : 11 ) ( xmin : -0.50 ) ( xmax : 10.50 )
( VALUE : 3 100 100 2 )
```

/ANALYS/DEFINITION/HST2

Define ZONE ,IDs of 2D-histograms to see. By defaults () press or value (0 given) the output is created with *single* histogram. See Section 2.3.2 for examples.

/ANALYS/DEFINITION/LIST1D

List 1D-histograms.

2.3. COMMANDS

```
ANAPAW> list1d [ HID ]
```

```
HID : 0, All lists are displayed.  
      -1, Definitions of histograms defined  
          by HST1 are displayed.
```

/ANALYS/DEFINITION/LIST2D

List 2D-histograms.

```
ANAPAW> list2d [ HID ]
```

```
HID : 0, All lists are displayed.  
      -1, Definitions of histograms defined  
          by HST2 are displayed.
```

/ANALYS/DEFINITION/LISTPROF

Print definition of *profile* histograms.

```
ANAPAW> listprof [ HID ]
```

```
HID : 0, All lists are displayed.
```

/ANALYS/DEFINITION/LISTGATE

Print definitions of the gates defined by ANA-CODE.

```
ANAPAW> listgate [ GKIND ]
```

```
GKIND : ' ' , All Gates are displayed  
        G [GATE] , Gates defined by [GATE] are displayed  
        A [AND] , Gates defined by [AND] are displayed  
        O [OR] , Gates defined by [OR] are displayed  
        X [XYGATE], Gates defined by [XYGATE] are displayed
```

/ANALYS/DEFINITION/AVIEW

Maybe the easiest way to view ANA-CODE definitions: gates ([GATE],[AND],[OR],[XYGATE]), and histograms (1D and 2D) see Section 3.4, is by opening a special X-Windows !

```
ANAPAW> aview [ OC ]
```

```
OC : 0 Open window (default)  
      C Close window
```


/ANALYS/DEFINITION/H CUT

Create 2D-gate data file. One needs to input gate ID before making gates. If $HID = 0$, current histogram is displayed. If filename=' ' the name **Anapaw.cut** is given. If you uses special options for the current histogram, see Section 2.3.3, you must add them here with histogram ID. Cut is made directly on window, using mouse ; left-click to add a new point and right-click to end the cut. Note that the last and the first point are automatically connected. Gates are mainly described in ANA-CODE, see Section 3.4.

```
ANAPAW> hcut [ HID filename COPT ]

          HID      : Histogram ID
          filename : File to save 2D-Gate Data
COPT      : M Main window
          S Scatter window
```

◇◇Example:

```
ANAPAW> hcut 117 bokuno.cut

ANAPAW-M : FileName = bokuno.cut

Input ID of 2D-Gate? (=1001)
Any Comment? > Wakatta !
```

Create a cut in the histogram number 117 and export it with the gate number 1001 in the **bunkuno.cut** file.

2.3.3 /ANALYS/HISTOGRAM/**/ANALYS/HISTO/HSTORE**

Store the Histograms in //PAWC to HBOOK type file. Default name file is **Anapaw.hbk** .

```
ANAPAW> hstore [ HBKFILENAME ]

          HBKFILENAME : Filename
```

◇◇Example:

```
ANAPAW> hstore text.hbk
```

/ANALYS/HISTO/FETCH

Fetch histogram ID from HBOOK file to //PAWC.

```
ANAPAW> fetch [ HBKFILENAME OFFSET ]

          HBKFILENAMEC : Filename to load the Histograms
          OFFSET      : Offset for Histogram ID
```

2.3. COMMANDS

ANALYS/HISTOGRAM/HLIST or I

Print list of all histograms.

```
ANAPAW> HLIST [ CHOPT ]
```

```
CHOPT : P Pager more
```

```
ANAPAW> I [ CHOPT ]
```

```
CHOPT : P Pager more
```

/ANALYS/HISTO/HT

Plot a histogram. This command saves current histogram ID. By default plot the current histogram. Options definition needs to put histogram ID.

```
ANAPAW> ht [ HID CHOPT ]
```

```
HID   : Histogram ID
```

```
CHOPT : Plot options, same as normal PAW
```

/ANALYS/HISTO/HTP

Plot the *current* histogram. Allow options definition without putting the histogram ID.

```
ANAPAW> htp [ CHOPT ]
```

```
CHOPT : Plot options, same as normal PAW
```

◆◆To plot the current histogram with colorscale enter:

```
ANAPAW> htp colz
```

/ANALYS/HISTO/HNT or /ANALYS/HISTO/HN

Plot the next histogram. This command saves current histogram ID.

```
ANAPAW> hn [ CHOPT ]
```

```
CHOPT : Plot options, same as normal PAW
```

/ANALYS/HISTO/HBT or /ANALYS/HISTO/HB

Plot the previous histogram. This command saves current histogram ID.

```
ANAPAW> hb [ CHOPT ]
```

```
CHOPT : Plot options, same as normal PAW
```

/ANALYS/HISTO/BLOW

Create a histogram between BLXMIN and BLXMAX.

```
ANAPAW> blow BLXMIN BLXMAX
```

```

          BLXMIN : Minimum
BLXMAX : Maximun
```

/ANALYS/HISTO/XBLOW

Expand a histogram using mouse in x-axis.

```
ANAPAW> xblow
```

/ANALYS/HISTO/XYBLOW

Expand a histogram using mouse in x-axis and y-axis.

```
ANAPAW> xyblow
```

/ANALYS/HISTO/PRX

Create a histogram for projection on x-axis. Note that this command does not only create a projection but does create the corresponding histogram.

```
ANAPAW> prx
```

/ANALYS/HISTO/PRY

Same as prx but for y-axis.

```
ANAPAW> pry
```

/ANALYS/HISTO/BNX

Create a histogram for projection on x-axis restricted to the y interval (YMIN, YMAX). Note that this command does not only create a projection but does create the corresponding histogram.

```
ANAPAW> bnx YMIN YMAX
```

```

          YMIN : Minimum
          YMAX : Maximun
```

/ANALYS/HISTO/BNY

Same as BNX but for the y-axis using restriction on x-axis.

```
ANAPAW> bnx XMIN XMAX
```

```

          XMIN : Minimum
          XMAX : Maximun
```

2.3. COMMANDS

/ANALYS/HISTO/SLX

Create projections onto the y axis, in x-slices. Note that this command does create the SNY corresponding histograms.

```
ANAPAW> slx SNY
```

```
SNY : Number of slices
```

/ANALYS/HISTO/SLY

Same a SLX but for the y-axis.

```
ANAPAW> sly SNX
```

```
SNX : Number of slices
```

/ANALYS/HISTO/MAMI

Change the y-axis limit of current histogram

```
ANAPAW> mami MRYMIN MRYMAX
```

```
MIN : Minimum
```

```
MAX : Maximun
```

/ANALYS/HISTO/XSTATUS

Print status of a given histogram

```
ANAPAW> xstatus
```

/ANALYS/HISTO/HSTATUS

Print status of the current histogram

```
ANAPAW> hstatus
```

/ANALYS/HISTO/ERASE

Erase histogram contents and **events summary** (EVS).

/ANALYS/HISTO/HDELETE

Delete the current (if HID1 = HID2 = 0) or given histograms.

```
ANAPAW> hdelete HID1 HID2
```

```
HID1 : First histogram
```

```
HID2 : Second histogram
```

/ANALYS/HISTO/XVAL; /ANALYS/HISTO/XYVAL

Print channel contents of a given bin in x-axis (XVAL) or in x and y axis (XYVAL).

```
ANAPAW> xyval [ CHOPT ]
```

```
CHOPT : ' ' Histogram mode
        S Scatter mode
```

/ANALYS/HISTO/HSTAT

Print statistics on *current* histogram.

```
ANAPAW> HSTAT
```

/ANALYS/HISTOGRAM/XSTAT

Print statistics on *selected* area for the *current* histogram.

```
ANAPAW> XSTAT
```

/ANALYS/HISTOGRAM/FIGA

Fit the *current* histogram using a Gaussian function. It then creates a new histogram.

```
ANAPAW> FIGA
```

/ANALYS/HISTOGRAM/CPRO

Created a new histogram using (contour) cut in another histogram, see Section [2.3.3](#).

```
ANAPAW> CPRO [ HID COPT ]
```

```
HID    : Histogram ID
CHOPT  : M Main Window
        S Scatter Window
        OPT HPLOT Option
```

/ANALYS/HISTOGRAM/AVY

Project contents to a profile histogram.

```
ANAPAW> AVY [ HID COPT ]
```

```
HID    : Histogram ID
CHOPT  : ' ' Default
        S Spread Option
        I Error on Mean
```

2.3. COMMANDS

/ANALYS/HISTOGRAM/XFITG

Fit the current histogram using a Gaussian function, but here you can select the abscissa area. It then creates a new histogram.

```
ANAPAW> XFITG
```

2.3.4 /GRAPHICS/

/GRAPHICS/PSON

Open a Postscript file for output.

```
ANAPAW> PSON
```

/GRAPHICS/PSOFF

Open the current output on Postscript file.

```
ANAPAW> PSOFF
```

◇◇Example with one histogram, sent to printer at the end.

```
ANAPAW> pson
File name (<CR>=)
-> ‘ ‘ sugoi.ps ’ ’
ANAPAW> ht 101
ANAPAW> psoff
ANAPAW> sh chsize sugoi.ps | lpr
```

Chapter 3

Analys Code

ANA-CODE aims to define gates and histograms for the carried experiment analysis. Everything is written in a file load in ANAPAW by the command book, see Section 2.3.2.

3.1 Hierarchy

Before starting with the code itself, one should know how the observables are organized in ANALYS. The first level is the analyzer type, which roughly correspond to the detector. Then the analyzer in divide in analyzer ID, like a detector can be divided in small parts. Finally these smaller parts are also divided in different words, which correspond to different observable (*raw or calculated*) like amplitude, timing, raw or calibrated data. . .




Figure 3.1: Analyzer's hierarchy in ANALYS

3.2 Comments

As Fortran file, the analysis description can be commented. The command looks like Fortran ones, that is to say using 'c'.

3.3. ANALYS

3.3 Analys

Analys File must start by the declaration of the analyzer used.

```
ANALYS
ID of the first analyzer
ID of the 2nd, etc
```

◇◇Example

```
ANALYS
2
73.4 Gate(s) 3.4.1 1D Histograms
```

3.4 Gates

3.4.1 GATE

For simple data gate, based on 1D-histogram:

```
GATE
gate ID, analyser type,
      ID of first analyzer into analyzer type
      ID of last analyzer into analyzer type
      word defining the analysis
      min X,max X
```

◇◇Example:

```
GATE
c TOF(RF2-F3) 54Ni
1, 2, 5, 5, 6, 263., 268.
```

3.4.2 Logical OR, AND

For logical gate OR (or AND):OR (or AND) gate ID, list of gates linked by OR (or AND), ...

◇◇Example:

```
AND
c 52Fe 850keV
200, 50, 903.5.
```

3.4.3 XYGATE

Create a gate using a cut file (ANAPAW cut file), see Section [2.3.2](#)

```
XYGATE
@,name of cut file
```

◇◇Example using a file:

```
XYGATE
c 1001
@,cut/54Fe.cut
```


3.4.4 STOP

Stop filling histogram placed after this command, for the current event.

```
STOP
[-]gate number
```

◇◇ Example with an already defined gate 3. In this case, all the commands placed after stop *won't* be read, if the current event does *not* satisfied (a minus is placed before) gate 3

```
GATE
3,1,4,4,5,25.,44.
STOP
-3
HST1
0,1,5,5,5,100,0,100.,'Beautiful 1D Histo gated'
```

3.5 Histograms

3.5.1 HST1: 1D Histograms

Histograms definition are equivalent with the HBOOK case. If you are not familiar with this CERN library please have a look at its manual:

<http://wwwasdoc.web.cern.ch/wwwasdoc/hbook.html3/hboomain.html>

The command to be put at the beginning of the 1D histogram list is HST1. Then each histogram is defined by:

```
HST1
gate ID, analyzer type,
      ID of first analyzer into analyzer in X
      ID of last analyzer into analyzer type
      word defining the channel read in X
      num of bin,min,max,'title'
```

3.5.2 PROFILE

This type of 1D histogram can be created using a 2D vector. The Y-coordinate is the *mean* of the vector's 2nd coordinate, with a error. The way to calculate this last depends of the last option, which can be ' ', 'S' or 'T'. Please see CERN HBOOK manual for further precisions.

Warning: the mean and the more probable 1st coordinate are most of time different. Please beware of your the 2nd coordinate distribution !

```
PROF
gate ID, analyzer type,
      ID of first analyzer into analyzer type in X
      ID of last analyzer into analyzer type in X
      word defining the channel read in X
      ID of first analyzer into analyzer type in Y
      ID of last analyzer into analyzer type in Y
```

3.6. ANA-CODE FILE EXAMPLE

```
word defining the channel read in Y
num of X bin,min X,max X,
min Y,max Y,
'title','error'
```

3.5.3 HST2: 2D Histograms

The command to be put at the beginning of the 2D histogram list is HST2. Then each histogram is defined by:

```
HST2
gate ID, analyzer type,
    ID of first analyzer into analyzer type in X
    ID of last analyzer into analyzer type in X
    word defining the channel read in X
    ID of first analyzer into analyzer type in Y
    ID of last analyzer into analyzer type in Y
    word defining the channel read in Y
    num of X bin,min X,max X,
    num of Y bin,min Y,max Y,
'title'
```

3.6 ANA-CODE file example

```
c === ANA-FILE ===
ANALYS
2
5
GATE
c TOF(RF2-F3) 54Ni oya
1, 2, 5, 5, 6, 263., 268.
c TOF(RF2-F3) 54Ni ko
2, 2, 5, 5, 6, 298.5, 303.5
c TOF(RF1-F3) 54Ni oya
3, 2, 5, 5, 4, 191.5, 196.5
c TOF(RF1-F3) 54Ni ko
4, 2, 5, 5, 4, 228., 232.
c TOF(RF2-F3) 52Fe oya
5, 2, 5, 5, 6, 270., 275.
AND
c 52Fe
11, 2, 4
OR
c SSD gate
21, 3, 4
c TOF
HST1
0, 2, 5, 5, 6, 300, 100., 400., 'TOF(RF2-F3)'
0, 2, 5, 5, 4, 300, 100., 400., 'TOF(RF1-F3)'
```

3.8. ANA-CODE FILE EXAMPLE29c SSD

3.6. ANA-CODE FILE EXAMPLE

```
HST2
c rot dE vs rotE
90,5,301,301,2,2,301,301,3,200, 1000.,4000.,200,-1000.,3000., 'dE vs E 2'
90,5,301,301,3,3,301,301,3,200, 1000.,4000.,200,-1000.,3000., 'dE vs E 3'
90,5,301,301,4,4,301,301,3,200, 1000.,4000.,200,-1000.,3000., 'dE vs E 4'
c
EXIT
```


Chapter 4

Anapaw Source Code

You can directly modify ANAPAW/Analys files, by modifying appropriate encoding files (`enc_*.f`). Each one describes the analyzer used (see Section 2.3.3) and can be modified, for example to calibrate the signal. The main used language is Fortran 77.

4.1 Analyzer Definition

To create a new analyzer, you should modify the source file named `usersrc.f` you can find in `usersrc-*.tar.gz`, see Section 2.1.2.

The minimum required is :

```
If(AnalyzerFlag(#))Then
  Call Enc_ANA( val(1,1,#),nx,ny,naok(#) )
EndIf
```

where `#` stands for the analyzer's number and `Enc_ANA` is a the encoding subroutine defined generally in a `enc_*.f` file.

`val(1,1,#)` represents the output array of observables values defined in the encoding subroutine, `nx` and `ny` give the two first dimension of `val(1,1,#)`.

Finally `naok` in an internal index that organize the observable in its array. Note that data are compressed: you do not need 100 lines in you data matrix `val` to access the word 100! The first column of you matrix contains this word ID....

An analyzer can be related to Raw Data. In this case, one should put in argument the corresponding array.

◇◇Example, where `EvtData(1)` is the raw data array, retrieved by the subroutine `FND_EID` and convert into observables in `Enc_SSD`:

```
c
c Analyzer 4: SSD
c
  If(AnalyzerFlag(4))Then
    call FND_EID(EvtData(1),5,iadr_sub_T,nw_sub_T) ! Timing
    call FND_EID(EvtData(1),6,iadr_sub_A,nw_sub_A) ! Analog
    If(nw_sub_T.le.0 .or. nw_sub_A.le.0)Then
```

4.2. ENCODING FILE

```
        EVTERR = .TRUE.
        Return
    EndIf
    Call Enc_SSD( EvtData(iadr_sub_T),nw_sub_T,
&              EvtData(iadr_sub_A),nw_sub_A,
&              val(1,1,4),nx,ny,naok(4) )
    EndIf
```

The purpose here is *not* to describe all the possible subroutines, but to give the main rules to modify them. Source files and their subroutines strongly depend on the experiment, the analysis, the author... Normally you will not have to hack this file, and if so, this manual may not enough to help you.

4.2 Encoding File

4.2.1 From raw data to observable

Here you can find the physics! Because you will certainly want to understand quickly how to deal with the code, we will here explain everything by examples.

Let's take 2 raw values given respectively by the index 7 & 11 of the array `raw_data`, i.e `raw_data(7)` & `raw_data(11)`. In 99% of the cases raw values are asked to be observable.

You will certainly see this type of code in the already written encoding files:


```
ID=2
naok=naok+1
val(1,naok) = ID
val(2,naok) = raw_data(7)
val(3,naok) = raw_data(11)
```

In this case ID is the IDentification for the sub-analyzer, as it is explained in Section 3.

Warning: You must write `val(1,naok) = ID`, because ANAPAW uses this first valu as identifier !

4.2.2 A first example

Let's take the case of long scintillator of length L with one photo tube at each side (let's say *left* and *right*), see Fig. 4.2.2.



Because the amplitude A received by the light is mainly exponentially decreasing with the distance x where a particle has hit the detector, the position

4.2. ENCODING FILE

of the hit can be guess by taking the logarithm of the ratio of the amplitude at each side. In equations:

$$A_{right} \propto e^{-\frac{L-x}{x_0}} \quad \& \quad A_{left} \propto e^{-\frac{x}{x_0}} \quad \Rightarrow \quad x \propto Ln \frac{A_{right}}{A_{left}}$$

In order to use the previous example, we can stand that the raw data for right and left amplitude are given by `raw_data(7)` & `raw_data(11)` respectively. We want this observable to be sort with the same sub analyser, in a encoding file named for example `enc_scinti`. Then the previous code becomes for example:

```
ID=2
naok=naok+1
val(1,naok) = ID
val(2,naok) = raw_data(7)
val(3,naok) = raw_data(11)
if (val(3,naok).ne.0) then
  pos_A = ln(val(2,naok)/val(3,naok))
  val(3,naok) = pos_A
endif
```

4.2.3 Calibration

Position of hits can also be known by time difference: $x \propto t_{left} - t_{right}$. Imagining that `raw_data(12)` & `raw_data(16)` are timing for right & left PMT, the previous code can be extended as:

```
ID=2
naok=naok+1
val(1,naok) = ID
val(2,naok) = raw_data(7)
val(3,naok) = raw_data(11)

c position using amplitude
if (val(3,naok).ne.0) then
  pos_A = ln(val(2,naok)/val(3,naok))
  val(3,naok) = pos_A
endif

val(4,naok) = raw_data(12)
val(5,naok) = raw_data(16)

c position using time
pos_t = val(5,naok)-val(4,naok)
val(6,naok) = pos_t
```

If we want to calibrate this values to get the position in mm, we need to hack 2 other files: one where the declaration of variables are (`*.inc`, here for example `scinti.inc`) and the source files for initializing the coefficients for calibration (`*_prm.f` here `scinti_prm.f`). Finally the coefficients themselves will be put in a `*.prm` file. The name of this file can be changed in the `$USER_SOURCE/setana` file, see Section 2.1.2.

4.2. ENCODING FILE

Starting from the end, we want to create a new sub analyzer in order to distinguish between raw and calibrated data.

```
ID=3
naok=naok+1
val(1,naok) = ID
val(2,naok) = to_x(1) + pos_A * to_x(2)
val(3,naok) = to_x(3) + pos_t * to_x(4)
```

where `to_x` is the array for calibration. We have to declare it in `scinti.inc` by adding a line to the common variables.

```
Real to_x
Common/scintillator/
.
.
.
& to_x(4)
```

In order to read the parameter file, for example, `scinti.prm`, you should first now under which number this last is open in `scinti_prm.f`. As an example let's take 80. The procedure to fill a array of real reading a file is:

```
READ_FLT(Unit_number,
         Array_name,
         Size,
         Error_code)
```

For an array of integer, the syntax is the same, but the function is :

```
READ_INT(Unit_number,
         Array_name,
         Size,
         Error_code)
```

The error code is different from 0 in case of problems ¹.

In this example, the code added is:

```
Call read_flt(80,to_x,4,ier)
```

Warning: The parameter are read continously in the `*.prm`, if you exchange 2 `read_flt` procedures without exchanging the parameters in the read file, then you will have no error message (in Fortran !) and the result will be...how shall I put it...stlightly weird.

As for the `scinti.prm` file:

```
c Calibration for scintillator
10, 0, .43, -7.4.3.
```

Of course a better way to hack this is to create two arrays: one for time, and one for analog. You can also create a condition for turning off the calibration. The final code will be :

```
◇◇scinti.inc
```

¹You can then add a condition to write a message, stop reading...


```

Integer iflag_AtoX_scinti,iflag_dttoX_scinti
Real A_to_x, dt_to_x
Common/scintillator/
.
.
.
& iflag_AtoX_scinti, A_to_x(2)
& iflag_dttoX_scinti, A_to_x(2)

◆◆scinti_prm.f

c Calibration for scintillator
  Call read_flt(80,A_to_x,4,ier)
  Call read_flt(80,dt_to_x,4,ier)

◆◆scinti.prm

c Calibration for scintillator : Analog
c iflag_AtoX_scinti
1
c values
10, 0
c Calibration for scintillator : Time
c iflag_dttoX_scinti
0
c values
.43, -7.

◆◆enc_scinti.f

ID=3
naok=naok+1
if (iflag_AtoX_scinti.ne.0) then
val(2,naok) = A_to_x(1) + pos_A * A_to_x(2)
endif
if (iflag_dttoX_scinti.ne.0) then
val(3,naok) = dt_to_x(1) + pos_t * dt_to_x(2)
endif

```

4.2.4 Multi array detector

Finally, in some cases, this plastic can belong to a bigger detector array. Imagining that our array is composed of 10 plastics, identicals, for which we want the same observable, then a quite look at how everything is organised in ANALYS tells us that we don't have to repeat each observable declaration, but loop on ID. and detector IDs will be taken from 1 to 10. Then the code becomes :

```

◆◆scinti.inc

Integer number_det
Parameter (number_det=10)
Integer iflag_AtoX_scinti,iflag_dttoX_scinti
Real A_to_x(2,number_det), dt_to_x(2,number_det)

```

4.3. NTUPLES

```
Common/scintillator/
.
.
.
& iflag_AtoX_scinti, A_to_x
& iflag_dttoX_scinti, A_to_x

◇◇scinti_prm.f

c Calibration for scintillator
  Call read_flt(80,A_to_x,2*number_det,ier)
  Call read_flt(80,dt_to_x,2*number_det,ier)

◇◇scinti.prm

c Calibration for scintillator : Analog
c iflag_AtoX_scinti
1
c values
10., 0., 26., 38., 18., 39., 38., 76., 43., 12.,
23., 62., 59., 53., 20., 23., 78., 23., 45., 89.
c Calibration for scintillator : Time
c iflag_dttoX_scinti
0
c values
.43 -7. .2 4.8 -45.6 -6.9 5.1 47. 90. 5.7
5.0 1.6 1.0 2.5 -2.5 -2.0 1.5 -2.0 3.4 2.8
c

◇◇enc_scinti.f

Do ID = 1,10
  naok=naok+1
  if (iflag_AtoX_scinti.ne.0) then
    val(2,naok) = A_to_x(1,ID) + pos_A(ID) * A_to_x(2,ID)
  endif
  if (iflag_dttoX_scinti.ne.0) then
    val(3,naok) = dt_to_x(1,ID) + pos_t(ID) * dt_to_x(2,ID)
  endif
Enddo
```

4.3 Ntuples

Anapaw allows the storage of ntuples (Column-Wise types), using the command ANTUPLE, see Section 2.3.1. But before using it, you should defined ntuples in two subroutines in `add_ntuple.f` file.

In subroutine NTBOOK you call the HBOOK routine HBNAME to describe the variables that are to be stored in the Ntuple. Please refer to HBOOK manual for more details, here under is the (simplified) synopsis:

```
CALL HBNAME (ID, CHBLOK, VARIABLE, CHFORM)
```

4.3. NTUPLES

ID
 Identifier of the Ntuple as in the call to HBNT.

CHBLOK
 Character variable of maximum length 8 characters
 specifying the name by which the block of variables
 described by CHFORM is identified.

VARIABLE
 The first variable that is described in CHFORM.

CHFORM
 Character string describing the variables to be
 stored in block CHBLOK

Warning: For now the ntuple number is automatically 10 and named EXP.
 In the second subroutine, the different variables will be filled, every other
 events.

◇◇ Here under a complete example:

```

c =====
c Subroutine NTBOOK
c -- from ANALYS Source !! Do Not Edit !! ----- c
c Logical AnalyzerFlag(50),INITENCFLAG(50)
c Logical USERFLAG(10),EVTERR
c Common/ANALYSLOGIC/ AnalyzerFlag,INITENCFLAG,USERFLAG,EVTERR
c ----- c
c -- for ntuple defined by user ----- c
c -- Analyzer 1
c Integer RunNum
c Integer Coin
c Common/CoinReg/ RunNum,Coin
c == Booking Part ==
c
c If(AnalyzerFlag(1)) Then
c Call HBNAME(10, 'CoinReg', RunNum,'RunNum:I,Coin:I')
c EndIf
c
c Return
c End
c =====
c Subroutine Add_Ntuple(IOFLAG)
c -- for Local ----- c
c Logical IOFLAG ! ( True : Next File )

```

4.3. NTUPLES

```
c -- from ANALYS Source !! Do Not Edit !! ----- c
  Real    val(500,500,50)
  Integer RunNumber
  Logical AnalyzerFlag(50),INITENCFLAG(50)
  Logical USERFLAG(10),EVTERR
  Common/ANALYSVALUE/  val
  Common/AnalysNtuple/ NTUPLE_FIRST,RunNumber
  Common/ANALYSLOGIC/  AnalyzerFlag,INITENCFLAG,USERFLAG,EVTERR
c ----- c

c -- for ntuple defined by user ----- c

c -- Analyzer 1
  Integer RunNum
  Integer Coin
  Common/CoinReg/ RunNum,Coin

c ----- c
  If(IOFLAG)Then                !! Do Not Edit !!
    Call Ntuple_io(0)           !! Do Not Edit !!
    Call Ntuple_io(1)           !! Do Not Edit !!
  EndIf                          !! Do Not Edit !!

c ----- c
c                               Definition of val(Word,ID,Analyzer)
c ----- c

c -- Analyzer 1
  If(AnalyzerFlag(1)) Then
    Coin = val(3,1,1)
  EndIf

  Return
  End
```

Chapter 5

Miscellaneous

5.1 Phone numbers

Room	Extension
Rips	4186
Control room	4126 or 4127
E6	4176
J5	4188

5.2 Data copy

5.2.1 Tape duplication

In order to create a backup, from one tape (in `/dev/nst0`) to another (in `/dev/nst1`)

```
mtcopy /dev/nst0 /dev/nst1
```

5.2.2 From tape to hard drive

In order to save data from tape (in `/dev/nst0`) on hard disk.

```
mtcopy /dev/nst0 /home/...
```

5.2.3 From hard drive to hard drive

Transferring data from the current hardware to, for example, `rips00` :

```
scp run[number].rdf rips00:[user]/[exp]
```

Reminder : most of time the current hard drive has a very limited size, so delete the data from it after transfer.

5.3. DATA STRUCTURE

5.3 Data structure

Block Header

```
0001 0000 0000 0000 Header Block
0000 0000 0000 0000 Data Block
```

Header Block 16KBytes

```
0001 0000 0000 0000 0000 0000 0000 --> Constant
0000 0000 ....
```

Word

```
0      : 0001 - Flag of Header
1-9    : 0000
10-13  : Run Number (ASCII 8Char.) "RUN-1066"
14     : Space * 2 (ASCII 2Char.) " "
15-23  : Start Time (ASCII 18Char.) "START => 12:58:56 "
24-32  : Stop Time (ASCII 18Char.) " STOP => 13:53:36 "
33     : Space * 2 (ASCII 2Char.) " "
34-42  : Print Time (ASCII 18Char.) "Print -> 13:57:07 "
43-47  : Print Date (ASCII 10Char.) " 22-SEP-98"
48-49  : 0000
50-89  : Header (ASCII 80char.)
90-129 : Space (ASCII 80Char.)
130-   : 0000
```

---- Example ----

```
\0 \0 \0 \0 R U N - 1 0 6 6 S T
A R T = > 1 2 : 5 8 : 5 6
S T 0 P = > 1 3 : 5 3 : 3
6 P r i n t = > 1 3 :
5 7 : 0 7 2 2 - S E P - 9 8
\0 \0 \0 \0 R u n 1 0 6 6 s a m e
a s b e f o r e
```

Ender Block 16KBytes

```
ffff 0000 0000 0000 0000 0000 0000 0000 -- > Constant
0000 0000 ....
```

Word

```
0      : ffff - Flag of Ender
1-9    : 0000
10-13  : Run Number (ASCII 10Char.) "RUN-1066"
14     : Space * 2 (ASCII 2Char.) " "
15-23  : Start Time (ASCII 18Char.) "START => 12:58:56 "
24-32  : Stop Time (ASCII 18Char.) " STOP => 13:53:36 "
33     : Space * 2 (ASCII 2Char.) " "
34-42  : Print Time (ASCII 18Char.) "Print -> 13:57:07 "
43-47  : Print Date (ASCII 10Char.) " 22-SEP-98"
48-49  : 0000
```

5.4. SCALER DATA

```

50-89 : Header      (ASCII 80char.)
90-129 : Ender      (ASCII 80Char.)
130-   : 0000

#### Event Block 16KBytes ####
0000 0000 0000 0000 .....

Word
0-3   : 0000 - Flag of Data
4     : Event size (Words - Include this)
5     : FID = 1 -> Constant
6     : Event ID
7     : Segment size (Words - Include this)
8     : Segment ID
9-    : Data

$ Event size $
Ex.) 0x8046 --> 1000|0000|0100|0110
           ^
           Always '1'   Event size

$ Fixed data (Non FERA) $
0       : Segment size (Words - Include this)
1       : Segment ID
3-     : Data

$ Variable data (FERA) $
0       : Segment size (Words - Include this)
1       : Segment ID
3       : FERA Data size (Words - Except this)
4       : FERA Header
Ex.) 0x800e --> 1000|1000|0000|1110
           ^
           Always '1'

16           11
-----
| 1|WC=0to15 |           |
-----

```

5.4 Scaler data

Scaler data is written by *babarDAQ* at the end of each block (16kB). The order is given by order of calling functions

```
read_scaler()
```

in *scaler.c* of the DAQ code. Each scaler value occupies 2 words (4 bytes) in block. For simplicity let's assume the following example: we have 12 scaler values to read in the DAQ code. Then the word order in block is the following:

Word#:	Data:
8168	Scaler1
8170	Scaler2
8172	Scaler3
....
....
8188	Scaler11
8190	Scaler12

There is 8192 words in a block so the word number of the first scaler value is always = $8192 - N * 2$. Where N is number of scaler values defined in the DAQ code.

Notes

